

	GÉNÉRALISTE			GRANDS ANIMAUX			ORGANES / MALADIES				
	Comprehensive Plus	Electrolyte	Pre-Surgical	Large Animal	Bovine	Equine	Kidney	Liver	Kidney Monitoring	Liver Monitoring	Diabetes
ALB	X		X	X		X	X	X			
ALP	X		X	X		X		X			
ALT	X		X					X		X	
AMY	X						X				
AST				X	X	X		X			
ALB/GLOB*	X		X	X		X		X			
BUN/CREA*	X		X	X		X	X		X		
BUN	X		X	X		X	X	X	X		
Ca	X			X	X	X	X				
CHOL	X										X
CK				X	X	X					
Cl <sup>-</sup>		X									
CREA	X		X	X		X	X		X		
FRU											X
GGT	X			X		X		X		X	
GLDH						X					
GLOB*	X		X	X		X		X			
GLU	X		X	X		X	X	X			X
K <sup>+</sup>		X									
LIP**	X										
Mg					X						
Na <sup>+</sup>		X									
NA <sup>+</sup> /K <sup>+</sup> *		X									
PHOS	X				X	X	X		X		
TBIL	X			X		X		X		X	
TP	X		X	X		X		X			
TRIG											X

\* Paramètres calculés.

\*\* Basé sur une excellente corrélation avec le réactif DGGR.